



Europäisches Patentamt
European Patent Office
Office européen des brevets



(11) **EP 1 426 008 A8**

(12) **CORRECTED EUROPEAN PATENT APPLICATION**

Note: Bibliography reflects the latest situation

(15) Correction information:
Corrected version no 1 (W1 A1)
INID code(s) 72

(51) Int Cl.7: **A61B 5/022**

(48) Corrigendum issued on:
06.10.2004 Bulletin 2004/41

(43) Date of publication:
09.06.2004 Bulletin 2004/24

(21) Application number: **03024942.9**

(22) Date of filing: **29.10.2003**

(84) Designated Contracting States:
AT BE BG CH CY CZ DE DK EE ES FI FR GB GR
HU IE IT LI LU MC NL PT RO SE SI SK TR
Designated Extension States:
AL LT LV MK

(72) Inventors:
• **Shirasaki, Osamu, Omron Healthcare Co., Ltd.**
Kyoto-shi Kyoto 615-00084 (JP)
• **Kubo, Takeshi, Omron Healthcare Co., Ltd.**
Kyoto-shi Kyoto 615-00084 (JP)

(30) Priority: **03.12.2002 JP 2002351136**


(74) Representative: **Kilian, Helmut, Dr.**
Wilhelms, Kilian & Partner
Patentanwälte
Eduard-Schmid-Strasse 2
81541 München (DE)

(71) Applicant: **Omron Healthcare Co., Ltd.**
Kyoto 615-0084 (JP)

(54) **Blood pressure measuring apparatus**

(57) A blood pressure measuring apparatus oppresses blood vessels of a human body by an inflating unit through a cuff. When the blood vessels are oppressed through the cuff, a pulse wave superposed on a cuff pressure signal through a cuff pressure detector is detected by a pulse wave detector. The waveform of the detected pulse wave changes similarly to changes of pressure in the oppressed blood vessels in one heart-beat period. The information of the detected pulse wave waveform and the cuff pressure at this time are stored in a memory. A blood pressure calculator determines the scale of pulse wave waveform on the basis of the information stored in the memory, and matches the determined scale of pulse wave waveform with the scale of waveform showing pressure changes in the blood vessels, so that a plurality of timings of coincidence of cuff pressure and blood pressure are detected, and the systolic pressure and diastolic pressure are calculated by using these detected information and stored information of pulse wave waveforms. It is therefore possible to measure blood pressures in a short time.

Fig.8



EP 1 426 008 A8